



Laser-based Headspace Inspection



Container Closure Integrity of Sterile Pharmaceutical Containers

Laser-based Headspace Analysis

Richard Millett
Vice President, Lighthouse Instruments Europe

Definition

Sterile product Container Closure Integrity (CCI)

The ability of a container (vial, ampoule, syringe, cartridge, bottle etc) to:

- Keep the contents **IN**
- Keep the contaminants **OUT**

Container Closure Integrity impacts.....



What is changing in the guidance & regulations?

- A new revised USP <1207> was implemented in August 2016
- EU Annex 1 undergoing revision
- FDA has announced a revision of their container closure guidance*

What does it mean for the Industry?

- Regulators increasingly critical of CCI data from legacy methods (blue dye, microbial ingress)
- Trend towards quantitative (deterministic) analytical methods
- Emphasis on Science-based justification
- Drive towards a coherent CCI strategy across the Product life cycle



USP <1207.1> Leak Detection Index

	Detectable Leaks	
Limit of Detection Index	Air Leak Rate (std cc/S)	Orifice Leak Size (um)
Class-1	$<1 \times 10^{-6}$	<0.1
Class-2	10^{-6} to 10^{-4}	0.1 to 1
Class-3	6×10^{-4} to 4×10^{-3}	2 to 5
Class-4	5×10^{-3} to 1.6×10^{-2}	6 to 10
Class-5	0.017 to 0.360	11 to 50
Class-6	>0.360	>50

Air leak rate at 1-atm differential pressure at 25 C, i.e. vial at full vacuum

USP <1207.2> Leak Detection Summary

	Leak Detection Class	Measurement Outcome
Tracer-gas	Class 1-4	Helium Loss
Laser-Headspace	Class 1-6	Gas Composition or Gas Pressure
HVLD	Class 3-6	Electrical Current
Pressure Decay	Class 3-6	Pressure Drop
Vacuum Decay	Class 3-6	Pressure Rise
Mass Extraction	Class 3-6	Mass Flow

While no single method is appropriate for all types of containers, Laser Headspace analysis is the only method for all sizes of defects

Characterizing the headspace non-destructively




What gases can be measured?

- Headspace oxygen
- Headspace carbon dioxide
- Headspace moisture (water vapor)
- Headspace total pressure levels





CCI testing – how?



Headspace Measurement

Exchange of gas between the container and the outside environment through a defect

CCI testing: other situations





What type of product-packages?

- Sterile liquid, or lyophilized, or dry-powder filled
- Transparent rigid containers:
 - Clear or amber glass
 - Transparent plastics
- Vials, syringes, ampoules, cartridges
- Nominal volume ranging from 0.2mL to 250mL



Advantages & disadvantages

Advantages


- Non-destructive
- Rapid
- Quantitative results
- Deterministic method
- Operator independent
- Applicable over whole leak range
- Permanent & temporary leaks detectable

Disadvantages

- Not all fill levels
- Sample needs to be transparent to laser
- Inline production inspection needs modified headspace



Equipment



100% inspection
Multiple heads for total headspace analysis.

Each system comes with calibration and reference standards prepared from your glassware.






Measurement performance

N=100	Headspace Oxygen (% atm)			
Standard Label	Known Value	Meas. Mean	Error	St. Dev.
0.0	0.000	0.01	0.01	0.02
1.0	1.005	0.96	-0.04	0.03
2.0	2.004	1.98	-0.03	0.03
4.0	3.998	4.02	0.02	0.04
8.0	7.999	8.13	0.13	0.03
20.0	20.00	19.93	-0.06	0.04

↑ Accuracy ↑ Precision





Part 2

GAS DIFFUSION THEORY



CCI testing: Gas diffusion theory

$$\vec{J} = -D\vec{\nabla}n \quad \text{Fick's 1st Law}$$

$$\frac{\partial P_i(t)}{\partial t} = \frac{-D \cdot A_0}{V} \frac{\partial P_i(z, t)}{\partial z}$$

New USP <1207> states:

“Mathematical models appropriate to leak flow dynamics may be used to predict the time required for detecting leaks of various sizes or rates.”

$$\%Oxygen(t) = 20.9\% \cdot \left[1 - \exp\left(-\frac{\alpha_{Diff}}{V} t\right) \right]$$

The change in oxygen concentration will be exponential with respect to time


Diffusion
Parameter

$$\alpha_{Diff} \left[\frac{cm^3}{s} \right] = \frac{D \cdot A_0}{L}$$

The Diffusion Parameter is a function of the Diffusion Coefficient, D , the defect cross-sectional Area, A_0 , and Depth, L .



Validation of Oxygen Ingress Model



With fixed values for:

$$D = 0.22 \text{ cm}^2/\text{s}$$

$$A_0 = 20 \mu\text{m}^2 \text{ (} 5 \mu\text{m } \varnothing \text{)}$$

$$V = 18 \text{cc (15R)}$$


Obtain an empirical
depth parameter value:

$$L = 6 \mu\text{m}$$

Model matches the data ± 0.3 %-atm oxygen at every point

Oxygen Ingress Model Example

Predicted oxygen concentration versus time for **ideal defects**





Part 3

INDUSTRY CASE STUDIES AND EXAMPLES



Overview of CCI case studies

1. Method development
2. Process optimization
3. Biologics Cold Storage CCI Study
4. 100% inspection of lyophilized product

Case study 1:


Method Development

Objective

- Detection of 5 micron leak within 30 minutes


Sample set

- 6R DIN clear tubing vial – 1.5mL product
- Positive controls: 2µm, 5µm, 10µm and 15µm laser drilled defects
 - Glass defects
 - Metal plate defects



Nominal hole size 5 µm

Image provided by Lenox Laser





Case study 1:

Method Development

– Phase 1: Manufacturing conditions

- Determine nitrogen purge conditions

– Phase 2: API reactivity

- Oxidation rate

– Phase 3: CCI Method development


- Diffusion test with vials with know defects (+ve controls)
- Effusion test with vials with know defects (+ve controls)
- Method protocol

Case study 1:

Method Development

Phase 1: Manufacturing conditions

- 50 product, water-filled and empty samples






Case study 1: Method Development

Phase 2: API reactivity

- 50 product samples opened to air and followed over time

**Mean measured headspace oxygen level
monitored over time**




	Oxygen (% atm)
Start	19.59
2 hours	19.50
26 hours	19.18
36 hours	17.63
144 hours	15.76

Case study 1: Method Development

Phase 3: CCI method development

- Diffusion tests with vials with known defects




Case study 1:

Method Development


Phase 3: CCI method development

- Effusion tests with vials with known defects



Case study 2:

Process optimisation



Case

- Liquid product in glass vial under N_2 atmosphere.
- All 200 vials passed visual inspection.

Result

- 192 accepted vials < 2% O_2
- 8 rejected vials \approx 20% O_2
- Total test time for 200 vials: < 45 minutes

Conclusion

Ineffective crimping caused defective vials with permanent leaks.



Case study 3: CCI testing for vials stored on dry ice (CO₂)

- 2R vials containing a Biologic, headspace 1 atm of air
- Stored on dry ice for 7 days.
- Thawed to room temperature (RT).
- Headspace conditions analyzed.
 - *Any change in the headspace condition would indicate a loss of CCI during deep cold storage*

Case study 3: CCI testing for vials stored on dry ice (CO_2)

- Air headspace vial at 1 atm at RT
- At low T, initial headspace condenses and creates underpressure
- Stopper can lose elastic properties & closure can be lost
- Cold dense gas from storage environment fills headspace
- Warming container to RT, stopper regains elasticity and reseals trapping the cold dense gas in the vial






Case study 3: CCI testing for vials stored on dry ice (CO₂)

- The cold dense gas trapped inside, expands as temperature increases creating overpressure
- Headspace gas composition could also change depending on storage environment
- Maintenance of changed headspace conditions can be monitored over time to verify that the leak was temporary.





Case study 3: CCI testing for vials stored on dry ice (CO₂)



Three different headspace measurements identify the same 3 vials as having CCI issues.




Case study 3: CCI testing for vials stored on dry ice (CO₂)

Some important comments on these results:

- Leaks during deep cold storage are usually temporary!
- CCI methods requiring an active leak (*blue dye, microbial ingress, pressure/vacuum decay*) will NOT identify these vials having temporary leaks.

Case study 4:

100% CCI testing of lyo product




Lyophilised
Product closed
at 200 mbar of
 N_2

Results of 6
consecutive
batches


Case study 4:

100% CCI testing of lyo product



- 1 atm air vials, gross (permanent) leaks
- Lyo headspace specified to be 200 mbar N₂
- If 800 mbar air enters vial = 16% O₂!
- Partial leaks stopped by capping

Case study 4: 100% CCI testing of Iyo product



Case 100% inspection

5 years of manufacturing data:

- 156 lots
- Total 1.9 million vials

Results

44-lots (28%) with zero rejects

Average reject rate was 0.25%

Difficult to manufacture a perfect CCI Iyo batch

Thank you for your attention



Demonstration

Let's consider the following product-package:

- Product ampoule closed at 500mbar N₂
- What are the headspace oxygen levels when the...



- A) ... container has retained CCl₄?**
- B) ... container has just lost closure?**
- C) ... container has permanently lost closure for some time?**

